Wide Product Range ... **PROMOTION** AUTUMN 2014

THE BRANDS YOU TRUST

Prices are recommended retail prices incl. VAT.

Promotions are valid until November 30, 2014 or as long as stocks last.

#COLTENE

SPECIAL RESTORATION

Hole in One

Filling in a single step with Fill-Up!™

Fill-Up!™

In a single step to give a perfect result

- Optimal depth polymerisation with minimal shrinkage due to dual curing system restoration
- Guaranteed single-layer technique even in very deep cavities of 10 mm
- Optimised sealing of margins reduced postoperative sensitivity
- Universal shade in a convenient Automix syringe for efficient placement

Deep. Fast. Perfect.

Fill-Up! INTRO KIT

6001 9342

 $2 \times 4.5 \,\mathrm{g}$ Fill-Up! Universal $1 \times 2.5 \,\mathrm{ml}$ Etchant Gel S

2 × 3 ml ParaBond Adhesive A & B

6 × Mixing Tip Short Super Fine

8 × Mixing Tip Short Fine

111 € + VAT

Fill-Up! ECONOMY KIT

6001 9355

5 × 4.5 g Fill-Up! Universal

 $15 \times Mixing Tip Short Super Fine$

20 × Mixing Tip Short Fine

199 € + VAT







Volume 26, No.2 SYNERGY D6

SYNERGY® D6 beauty made simple

- Ideal handling
- High aesthetics
- Excellent polishability

SYNERGY D6 SYRINGES + ONE COAT 7.0 WELCOME KIT

3 Syringes (A1/B1, A2/B2, A3/D3), Shade Guide, 5 ml One Coat 7.0, 50 microbrushes, accessories + SwissFlex Trail kit sample

99,80 € + VAT

3244

3243

6001 4361

6001 4362

SYNERGY D6 TIPS + ONE COAT 7.0 WELCOME KIT

20 Tips, (5 per A1/B1, A2/B2, A3/D3, Enamel Universal), Shade Guide, 5 ml One Coat 7.0, 50 microbrushes, accessories + SwisFlex Trail kit sample

99,80 € + VAT

SYNERGY D6 LIMITED EDITION 13 SYRINGES

5 Syringes each (A2/B2, A3/D3) 1 Syringes (A3.5/B3, A4/C3, Enamel Universal), Shade Guide

249 € + VAT

SYNERGY D6 LIMITED EDITION 130 TIPS

50 Tips each (A2/B2, A3/D3) 10 Tips each (A3.5/B3, A4/C3, Enamel Universal), Shade Guide

249 € + VAT

Dr. med. dent. Mario J. Besek Swiss Dental Center, Zürich, Switzerland

In contrast to commercially available "aesthetic composites", SYNERGY D6 is a simplified system for routine use that requires less effort, yet produces comparable or better results. Six dentine shades are available, which blend together so that fewer individual shades are required. The secret is in the chameleon effect with a balanced opacity and homogenised translucency, which are used for aesthetically reproducing the shade of the natural tooth shade. Two enamel shades – universal and opalescent – complete the system and can also be layered as required.

The working consistency has been optimally selected for the operator in all the required parameters that are important during sculpting and finishina

In summary, SYNERGY D6 is a viable, all-purpose practice system that successfully produces aesthetic results in the anterior and posterior region using a simple technique.



COMPONEER®

Surprise your patients with a new smile – in only one session!

- No laboratory required
- One session
- Naturally aesthetic corrections using freehand technique
- Easy application with prefabricated composite veneers
- Brilliant result
- Attractive added value

COMPONEER Intro System-kit

TIPS 6001 4221 SYRINGES 6001 4222

24 COMPONEER:

every OK Set: 11, 12, 13, 21, 22, 23

1 COMPONEER Set OK M – Enamel Universal

- 1 COMPONEER Set OK L Enamel Universal
- 1 COMPONEER Set OK M Enamel White Opalescent
- 1 COMPONEER Set OK L Enamel White Opalescent
- 1 COMPONEER accessories / Material
- 1 COMPONEER accessories / Instruments

893 € + VAT









For information on COMPONEER workshops please see www.componeer.info.com

COMPONEER® REFILL

COMPONEER Sets und Refills

Enamel White Enamel Universal Opalescent

COMPONEER Set **OK S, M, L, XL** 6 Veneers: 11, 12, 13, 21, 22, 23 COMPONEER Refill **OK S, M, L, XL** 2 Veneers

COMPONEER Set **OK Premular S, L** 2 Veneers

COMPONEER Set **UK S, M** 6 Veneers: 31, 32, 33, 41, 42, 43 COMPONEER Refill **UK S, M** 2 Veneers



Clinical Indications

Treatment of caries

- Extensive cervical and interdental caries, discolouration of enamel, and slight malposition
- 2. Caries treated and homogenised anterior with COMPONEER





Optimisation of old restorations

- 1. Extensive old and insufficient restorations
- 2. Reconstructed function with COMPONEER





Attrition, abrasion, erosion

- 1. Edge-to-edge occlusion with abraded and broken anterior teeth with old restoration
- 2. Functional and aesthetic teeth with normal overjet constructed with COMPONEER





Tooth fracture

- 1. Loss of 1/3 of anterior tooth by fracture
- 2. Simplified shaping and construction with COMPONEER





Aesthetic / cosmetic Indications

Extending incisors

- 1. Slightly shortened incisal edge with moderate malposition
- 2. The tooth proportions are reconstructed and aligned with COMPONEER





Malposition

- 1. Extensive malposition waith crowding
- 2. Axially aligned correction and repositioning with COMPONEER







SPECIAL ENDODONTICS

HYFLEX® CM

6001 1048 - 6001 1079

The new regenerative NiTi File

- No shape memory + Extreme flexibility = Superior canal tracking
- Regains shape after heat treatment = Multi-use
- Up to 300 % more resistance to separation

64,60 € + VAT



LIKE A PHOENIX FROM THE ASHES...
...THE REGENERATIVE NITI FILE!

No shape memory

- + Extreme flexibility
- = Superior canal tracking

HyFlex™ Controlled Memory NiTi Files have been manufactured utilizing a unique process that controls the material's, making the files extremely flexible but without the shape memory of conventional memory NiTi Files. This gives the file the ability to follow the anatomy of the canal very closely, reducing the risk of ledging, transportation or perforation. Further-more the files can be pre-bent, similar to the way of stainless steel. Particularly in root canals with abrupt curvatures this can help to avoid the creation of steps.

Controlled Memory








Up to 300 % more resistance to separation

Cyclic Fatigue Test*



CANALPRO CL

Cordless Handpiece

- Thin ergonomic design for better access to hard to reach areas
- Large, user-friendly LCD panel
- Memory for up to 5 programs
- Auto-reverse modus
- 6 head positions allows of contra angle head to be adjusted to your preferred position
- Operates continuously for approx. 80 minutes at rated load
- Quick charging time of approx. 90 minutes non-contact charger prevents improper charging due to deterioration of metallic terminal
- Energy-saving, turns off when not used for 10 minutes
- Contra angle head is autoclavable at 135°C
- Rotation Speed (min-1): 140 ~ 550
- Torque Level (N·cm): 0,3 ~ 3,0







6001 4334

Innovative and accurate

- Accuracy of measurement
- User-friendly, intelligent interface
- Virtual Apex Function
- Built-in self-test for a fast function control

1.014 € + VAT







SPECIAL PROSTHETICS

AFFINIS®

Perfect Impressions

- Surface affinity and optimal flow properties
- Minimizes the risk of air entrapments and marginal distortions
- Excellent readability

20 %

AFFINIS System 360

AFFINIS System 360 MonoBody Refill	6482
AFFINIS System 360 heavy body Refill	6487
AFFINIS System 360 putty Refill	6472

AFFINIS putty soft/super soft

AFFINIS putty soft Single Pack	6530
AFFINIS putty super soft Single Pack	6535

AFFINIS PRECIOUS System 50

AFFINIS PRECIOUS light body 50 ml	6775
AFFINIS PRECIOUS regular body 50 ml	6776

AFFINIS System 50

AFFINIS light body 50 ml	6501
AFFINIS regular body 50 ml	6511

COMPRECORD

- Air texturised retraction cord
- Outstanding absorbency





SPECIAL TREATMENT AUXILIARIES

HYGENIC® COMPLETE KIT

A dry path to success

H02778

HYGENIC Dental Dam is made of pure, natural rubber latex. Tough and hard to tear, the powdered material is used as a barrier during operative and endodontic procedures.

36 sheets:

Dental Dam (green, medium)

Plastic Frame

Fiesta color coded clamp pack with 9 clamps

Dental Dam Punch

Dental Dam Forceps

Template

Wedjets small

389 € + VAT



HYGENIC® ELASTI-DAM

Elasti-Dam is considerably more elastic1) and tear resistant than standard latex dam making application much easier. At the same time Elasti-Dam offers the excellent retraction properties valued in standard latex dam. Elasti-Dam is powderfree2) and has a controlled low level of protein, so patient and dental team are less exposed to allergenic particles and face a lower risk of developing a latex allergy.

H11625 Elasti Dam Medium green
H11627 Elasti Dam Medium blue
H11676 Elasti Dam Heavy green
H11690 Elasti Dam Heavy blue

20,30 € + VAT



HYGENIC® WEDJETS Stabilizing Cord

Alternative or in addition to conventional clamp method of securing dental dam.

Latex available in 3 sizes and non-latex in 2 sizes.

H08040 Latex Extra small (2,13 m) **H06522** Latex small (2,13 m)

H06523 Latex large (2,13 m)

H06524 Non-Latex small (2,13 m) **H06525** Non-Latex Large (2,13 m)

20 € + VAT



SPECIAL ROTARY INSTRUMENTS

DIATECH® SPEEDSTER

DIATECH Speedster carbide burs have an advanced technology design for rapid reduction of all materials (amalgam, gold, porcelain, enamel, PFM = Porcelain Fused to Metal). The computer designed flutes cut more concentric, eliminate chatter and vibration, cut longer and let you reduce the preparation and chair time!

Proprietary end cut design is excellent for Endo access

- Efficacité de coupe supérieure
- Durée de vie plus longue
- Carbure de tungstène monobloc concentrique



Which shape for which indication?





S1

Round

Used for initial cavity preparations, retention grooves, and to create an initial depth mark for prosthetic preparations.





S2

Inverted Cone / Double Cone

The shape of this carbide makes it ideal for removal of amalgams. It may also be used to contour the occlusal anatomy of final restorations.



S3

Pear

Recommended to contour the occlusal anatomy, cavity preparation and amalgam removal.





CA

Cross Cut Fissure

This carbide may be used to create butt-joint margins and to enhance the proximal box of preparations for ceramic inlay/onlays and crowns.





S5

Round End Cross Cut Fissure

Ideal for crown removal, creating deep chamfer crown preparations for ceramic restorations. Larger diameters may be used to create occlusal depth grooves.







S6

Cross Cut Tapered Fissure

Recommended for use in preparing shoulder margins and inlay/onlay preparations. May also be used to refine the axial walls of a preparation and to bevel or contra-bevel finish lines.



PROMOTION PROMORDIANT PROMOTION PROM

Coltène/Whaledent AG
Feldwiesenstr. 20
9450 Altstätten/Switzerland
Tel. + 41 71 757 53 00
Fax + 41 71 757 53 01
info.ch@coltene.com

For detailed information please contact:

Marzena Kotynska | +48 607 989 535

Marzena.Kotynska@coltene.com

Agnieszka Witamborska
Agnieszka.Witamborska@coltene.com

Slaskie-swietokrzyskie- podkarpackiemałopolskie | Tel. 0 795 159 515

KOL-Dental | Tel. 0 505 242 758

POLDENT Sp. z o.o. | Tel. (22) 351 76 50 do 55

